

Pipette Clinic

Homer Tsai


Overview

1. Principles in liquid handling
2. Guidelines to proper handling of pipettes
3. Decontamination and cleaning of pipettes

eppendorf
In touch with life

2

Liquids category



eppendorf
In touch with life

3

Principles in liquid handling

Pipetting principles

Air-cushion
(air displacement)

- *Forward pipetting*
- *Reverse pipetting*


Positive displacement



eppendorf
In touch with life

4


Air-cushion pipette



No contact
between
piston and sample


eppendorf
In touch with life

Air-cushion principle




eppendorf
In touch with life

Air-cushion pipettes - Manual




eppendorf
In touch with life

Air-cushion pipettes - Electronic



eppendorf
In touch with life

Forward pipetting – Technique



Press to **1st stop**.
Immerse tip a few
mm into liquid.

Release button
slowly. Tip fills up.

Press to **2nd stop**
to dispense.

eppendorf
In touch with life


Forward pipetting - Applications

Application is recommended for standard *aqueous* solutions:

- ❖ Water
- ❖ Buffer
- ❖ Diluted saline solutions
- ❖ Diluted acids and alkalis

eppendorf
In touch with life

Reverse pipetting – Technique



Press to **2nd stop**.
Immerse tip 3-4 mm
into liquid.

Release button
slowly. Tip fills up.

Press to **1st stop** to
dispense. Some liquid
will remain in the tip.

eppendorf
In touch with life


Reverse pipetting – Applications

Application is recommended for:

- ❖ High vapor pressure liquids
- ❖ Viscous liquids

eppendorf
In touch with life

High vapor pressure liquids




Viscous liquids

- Problem:
- High resistance to flow
 - Thin film of liquid remains in the tip

SOLUTION

- ✓ Slow aspiration and dispensing speed
- ✓ Reverse pipetting
- ✓ Pre-wetting

High density liquids



SOLUTION


- ✓ Re-adjustment of air-cushion pipette →
- Attention: adjustable-volume pipette → fix-volume pipette

Pipetting detergent-containing liquids

Following dispensing detergent-containing liquids, a film of liquid remains on the tip surface.

- Loss of sample
- Increased reagent consumption, e.g. qPCR master mix
- Poor precision (reproducibility)

- 1 Substantial sample loss with standard tips
- 2 Maximum sample recovery with **epTIPS LoRetention**



epT.I.P.S LoRetention

- The surface is rendered ultra-hydrophobic by modification at the molecular level without any
 - coating or
 - additives
 which may contaminate the sample
- Minimizes loss of samples in solutions containing detergents (e.g. Triton X-100, SDS, Tween)
 - ✓ maximum sample recovery
 - ✓ improved reproducibility

Tips to further improve pipetting accuracy

Holding angle

90°

90°

30° - 40°

Immersion depth

1 cm

3 cm

3-4 cm

Inaccuracy

0.2 - 0.4%

0.6 - 0.8%


1 - 1.2%

Principles in liquid handling

Pipetting principles

- Air-cushion (air displacement)
 - Positive displacement
- Forward pipetting
 - Reverse pipetting

Positive displacement principle




Positive displacement – Applications

- ✓ Ideal for “problematic liquids”:
 - High density
 - High vapor pressure
 - Viscous
 - ⇒ High accuracy and precision
- ✓ Contamination free – prevention of aerosols:
 - PCR
 - Toxic, infectious or radioactive liquids
 - Aggressive liquids (e.g. concentrated acids)

eppendorf
In touch with life

Pipetting viscous liquid

Effect of viscous liquid on pipetting accuracy and precision




21

22

Positive displacement dispensers

Multipette stream & Xstream


Pipetting
Dispensing
Automatic dispensing
Sequential dispensing
Aspiration
Titration




eppendorf
In touch with life

Contamination-free pipetting

Air-cushion pipette



Positive displacement with integrated piston in the tip



eppendorf
In touch with life

23

24

Minimize aerosol contamination with filter tips



Blue layer binds smaller aerosols (e.g. biomolecules)
Retains drops, splashes and larger aerosols


Double protection !



Hydrophobic surface

eppendorf
In touch with life

Filter tips efficiency – Smoke test



* See "Smoke Test" video

eppendorf
In touch with life

Handling problematic liquids

Liquids category	Air-cushion pipette			Positive displacement pipette
	Reverse pipetting	Pre-wetting	Slow aspiration and dispensing speed	
High vapor pressure	✓	✓		✓
Viscous	✓	✓	✓	✓
High density	Re-adjustment of pipette			✓
Detergent-containing	epTIPS LoRetention			
PCR/ Infectious/ Radioactive	Dualfilter tips			✓

eppendorf
In touch with life

Overview

- Principles in dispensing of liquids
- Guidelines to proper handling of pipettes

eppendorf
In touch with life

Guideline #1

**Know the liquid's physical properties,
choose the best pipetting techniques**
(forward, reverse or pre-wetting of tip)

- 1. **Increases** accuracy and reproducibility of results
- 2. **Reduces** contamination risks



Guideline #2

Never lay a pipette down with fluid in its tip!

Prevent fluid from entering the pipette
internal surfaces through the nose cone



Reduce risk of cross-contamination

Hold a pipette with a filled tip upright



Guideline #3

DO NOT place pipettes onto the surface of the bench

- Decrease contamination risks from chemicals that have been spilled onto the bench.
- pipette cone may get chipped/bent out of shape, resulting in poor tip fitting.

Use...



Pipette
carousel



Wall
mount



Guideline #4

Never release the push-button rapidly during aspiration



High aspiration speed promotes aerosol
formation and leads to contamination.

Use filter tips – forms a barrier between the
sample and the interior of the pipette.



Guideline #5

Accidentally absorbed liquids
SHOULD NOT be left to dry



Immediately clean the pipette

33

eppendorf
In touch with life

Guideline #6

To ensure continued precision and accuracy,
send the pipette for servicing and calibration.

The frequency of servicing and calibration depends on:

1. Frequency of pipette usage
2. Number of users

...but at least **once a year**



34

eppendorf
In touch with life

Overview

1. Principles in dispensing of liquids
2. Guidelines to proper handling of pipettes
3. Decontamination and cleaning of pipettes

35

eppendorf
In touch with life

Contamination scenarios


1. Aqueous solutions, buffers, acids and alkalis
2. Organic solvents
3. Potentially infectious liquids and cell cultures
4. Nucleic acids
5. Proteins
6. Radioactive substances

36

eppendorf
In touch with life

Contamination scenario 1

Aqueous solutions, buffers, acids and alkalis




Contamination scenario 2

Organic solvents

Disassemble pipette lower parts → Allow liquid to evaporate

Alternatively:



• DO NOT use concentrated detergent.

Contamination scenario 3

Potentially infectious liquids and cell cultures

Method A:

alcoholic disinfectants



Method B:

ultra-violet (UV) irradiation

Method C:

autoclaving



Contamination scenario 4


Nucleic acids

Method A:

Sodium hypochlorite (5-10%)

Method B:

Glycine/HCl buffer (pH 2.0)



Sodium hypochlorite (5-10%)

Use 5-10% fresh solution
↓
Soak contaminated lower parts
for 20 – 30 min
↓
Rinse well with dH₂O
↓
Allow to dry
↓
Lubricate piston
↓
Assemble pipette




41

eppendorf
In touch with life

1X Glycine/HCl (pH 2.0)

Boil contaminated lower parts
for ~ 30 min
↓
Rinse well with dH₂O
↓
Allow to dry
↓
Lubricate piston
↓
Assemble pipette




42

eppendorf
In touch with life

Contamination scenario 5

Proteins

Rinse contaminated lower parts
with mild detergent
↓
Rinse well with dH₂O
↓
Allow to dry
↓
Lubricate piston
↓
Assemble pipette



DO NOT use alcohol → insolubility of proteins

43

eppendorf
In touch with life

Contamination scenario 6

Radioactive substances

Rinse contaminated lower parts with
radioactive cleaning solutions
↓
Rinse well with dH₂O
↓
Allow to dry
↓
Lubricate piston
↓
Assemble pipette



44

eppendorf
In touch with life



Thank you !

Email: homertsai@hi-point.com.tw

eppendorf
In touch with life